



# National Institute of Standards & Technology

## Report of Investigation

### Reference Material 8321

#### Peptide Mixture for Proteomics

This Reference Material (RM) is intended to support investigations used to identify peptides in complex peptide mixtures such as those in mass spectrometry-based proteomics. RM 8321 can be used to help assess the confidence of peptide identification within a laboratory or comparability between laboratories or among different measurement approaches. RM 8321 can also be used in the development and validation of new investigative approaches for identifying peptides in complex peptide mixtures. A unit of RM 8321 consists of three vials, each containing approximately 50  $\mu\text{L}$  of frozen aqueous solution containing 0.1 mL/L formic acid. The peptides in this RM are estimated to be in a concentration range of 0.1 pmol/ $\mu\text{L}$  to 10 pmol/ $\mu\text{L}$ .

RM 8321 is an aqueous solution of approximately 440 synthetic peptides, present at a range of concentrations that span approximately three orders of magnitude. RM 8321 was designed to provide a complex mixture of peptides for evaluating the performance of proteomics mass spectrometry instruments coupled to liquid chromatography (LC) [1]. Peptides were chosen to cover the chromatographic “space” of a typical reverse phase gradient elution analysis, offering a range of elution profiles. The synthetic peptides in RM 8321 have the same amino acid sequence as tryptic peptides from 50 high abundance human plasma proteins that have been observed as proteotypic peptides through multiple published investigations by the proteomics community. Proteotypic peptides are those peptides which are observed repeatedly by mass spectrometric-based proteomics investigations by different investigators. Therefore, proteotypic peptides are expected to be readily released from enzymatic digestions of the precursor protein, ionize well by electrospray ionization, generate high quality tandem mass spectrometry (MS/MS) spectra, and are stable during the processes of sample preparation and analysis.

**Reference Values:** A NIST reference value represents the best estimation of the true value based upon the available data [2]. Table 1 lists a set of heuristic rules which describe the confidence of peptide identification in RM 8321. Table 2 lists the peptides present, grouped by confidence level, as determined by two types of LC-MS/MS analyses and comparison to mass spectral libraries [3].

**Expiration of Reference Values:** RM 8321 is valid, within the specified confidence levels, until **05 June 2020**, provided the RM is handled and stored in accordance with the instructions given in this report (see “Instructions for Storage and Use”). This report is nullified if the RM is damaged, contaminated, or otherwise modified.

**Maintenance of RM:** NIST will monitor this RM over the period of its validity. If substantive technical changes occur that affect the value assignment before the expiration of this report, NIST will notify the purchaser. Registration (see attached sheet or register online) will facilitate notification.

Overall direction and coordination of technical measurements leading to peptide identification were performed by D.M. Bunk of the NIST Biomolecular Measurement Division.

Analyses were performed by A.S. Beasley, M. Lowenthal, and D.M. Bunk of the NIST Biomolecular Measurement Division.

Support aspects involved in the issuance of this RM were coordinated through the NIST Office of Reference Materials.

Michael J. Tarlov, Chief  
Biomolecular Measurement Division

Gaithersburg, MD 20899  
Report Issue Date: 28 December 2015

Robert L. Watters, Jr., Director  
Office of Reference Materials

**Reference Values:** Table 1 defines the corresponding confidence levels of peptide identification in RM 8321. Table 2 identifies the peptide content, grouped by confidence level, based on both the data dependent acquisition (DDA) and dynamic multiple reaction monitoring (dMRM) analyses.

Table 1. Definitions of Heuristic Rules

| <i>Confidence Level</i> | <i>Heuristic Definition</i>  |
|-------------------------|--|
| High Confidence         | Peptides listed were observed in more than 90 % of all laboratory investigations (including technical and process replicate analyses). |
| Confident               | Peptides listed were observed in more than 50 % and less than 90 % of all laboratory investigations.                                   |
| Low Confidence          | Peptides listed were observed in less than 50 % of all laboratory investigations.  |

Table 2. Reference Peptide Content of RM 8321

**High Confidence:**

|                    |                    |                    |
|--------------------|--------------------|--------------------|
| AATVGSLAGQPLQER    | HLSLLTTLNLR        | SDVVYTDWKKDK       |
| ADLFYDVEALDLESPK   | HLVPGAPFLQALVR     | SEETKENEGFTVTAEGK  |
| AEDHFSVIDFNQNR     | HPDYSVLLLLR        | SELEEQLTPVAEETR    |
| AEFAEVSK           | HQLYIDETVNSNIPTNLR | SELTQQLNALFQDK     |
| AFQPFFVELTMPYSVIR  | HQTVPQNTGGKNPDPWAK | SFFSFLGEAFDGAR     |
| AGALNSNDAFVLK      | HSIFTPETNPR        | SGAQATWTELPWPHEK   |
| AGDFLEANYMNLQR     | HSTIFENLANK        | SGFPQVSMFFTHTFPK   |
| AGKEPGLQIWR        | HTLNQIDEVK         | SGKDPNHFRPAGLPEK   |
| AHVSFKPTVAQQR      | HTSVQTTSSGSGPFTDVR | SGVQQLIQYYQDQK     |
| AHYGGFTVQNEANK     | HVVPNEVVVQR        | SHALQLNLR          |
| AIEDYINEFSVR       | HYEGSTVPEK         | SIEVFGQFNGK        |
| AIGYLNTGYQR        | HYQINQQWER         | SKEFQLFSSPHGK      |
| AKPAEDLR           | IADAHLDR           | SKEQLTPLIK         |
| ALDLINKR           | IADNKQSSFK         | SLAELGGHLDQQVEEFR  |
| ALLVGEHLNIIVTPK    | IAFSATR            | SLAPYAQDTQEK       |
| ALMDETMK           | IAQWQSFQLEGGLK     | SLHTLFGDK          |
| ALTDMPQMR          | IDTQDIEASHYR       | SPELQAEAK          |
| ALVQQMEQLR         | IEGNLIFDPNNYLPK    | SPELQAEAKSYFEK     |
| ALYLQYTDETFR       | IHWESASLLR         | SSALDMENFR         |
| AMAVEDIISR         | IIRSSDPNEDIVER     | SSEDPNEDIVER       |
| APNHAVVTR          | IIVPLNLR           | SSLSVPYVIVPLK      |
| AQLVDMK            | IKVLNQELR          | SSNLIILEHLK        |
| AQRQVVAGLNR        | ILGGHLDK           | SVLGQLGITK         |
| ASEAEDASLLSFMQGYMK | IPIEDGSGEVVLSR     | SVNDLYIQK          |
| ASSIIDELFQDR       | IPLDLVPK           | SVSDGIAALDLNAVANK  |
| ASSIIDELFQDRFFTR   | ISASAEELR          | SVVDENFSWYLEDNIK   |
| ASTPNGYDNGIHWATWK  | ISEGLPALEFPNEK     | SYFEKSKEQLTPLIK    |
| ATEHLSTLSEK        | ITENDIQIALDDAK     | SYFPESWLWEVHLVPR   |
| ATFQTPDFIVPLTDLR   | ITPNLAEFASFSLYR    | SYTITGLQPGTDYK     |
| ATGVLYDYVNK        | IVSSAMEPDR         | TAAQNLYEK          |
| ATVVYQGER          | IYGNQDTSSQLKK      | TAAQNLYEKTYLPAVDEK |
| AVMDDFAAFVEK       | IYHSHIDAPK         | TAGWNIPMGLLYNK     |
| AVSMPSFSILGSDVR    | IYISGMAPRPSLAK     | TDAPDLPEENQAR      |
| AYKSELEEQLTPVAEETR | IYLYTLNDNAR        | TEDTIFLR           |
| DALSSVQESQVAQQR    | KAMAVEDIISR        | TEGDGVYTLNDK       |
| DAQYAPGYDKVK       | KATVVYQGER         | TEHPFTVEEFVLPK     |
| DDEEFIESNK         | KDNEQHVFVK         | TEHYEEQIEAFK       |
| DDNPNLPR           | KELSSFIDK          | TELRPGETLNVNLLR    |
| DFHINLFQVLPWLK     | KFPSGTFEQVSQVLK    | TEVNVLPQAK         |
| DFVQPPTK           | KGWVALNPLR         | TGAQELLR           |
| DGAGDVAFVK         | KLSSWVLLMK         | TGLQEVEVK          |
| DGNTLTYYR          | KLVPFATELHER       | THLAPYSDELK        |
| DHAVDLIQK          | KLWAYLTINQLLAER    | THLPEVFLSK         |

**High Confidence (continued):**

|                    |                    |                     |
|--------------------|--------------------|---------------------|
| DINYVNPVIK         | KSASDLTWDNLK       | THLPEVFLSKVLEPTLK   |
| DISEVVTTPR         | KTLLSNLEEAKK       | TILGTMPAFEVSLQALQK  |
| DKVNSFFSTFK        | KVEQAVETEPEPELR    | TLEAQLTPR           |
| DLATVYVDVLK        | KVPQVSTPTLVEVSR    | TLLPVSKPEIR         |
| DLKVEDIPLAR        | KWQEEMELYR         | TLLSNLEEAK          |
| DLMEKVKSPELQAEAK   | KYFIDFVAR          | TMEQFTIHLTVNPQSK    |
| DMYSFLEDMGLK       | KYNSQNQSNNQFVLYR   | TMTIHNGMFFSTYDR     |
| DRVVEESELAR        | LANLTQGEDQYYLR     | TNFDNDIALVR         |
| DSGRDYVSQFEGSALGK  | LAVYQAGAR          | TPSAAYLWVG TGASEAEK |
| DSQEEKTEALTSK      | LDAQASFLPK         | TSNFNAAISLK         |
| DSQEEKTEALTSKR     | LDEL RDEGK         | TTIEKPVWLGLGPIIK    |
| DTVIKPLLVEPEGLEK   | LDEVKEQVAEVR       | TTNIQGINLLFSSR      |
| DTVQIHDITGK        | LDGKFSVVYAK        | TVGSDTFYSFK         |
| DVVLFEK            | LDGSVDFK           | TVIGPDGHK           |
| DWHGVPGQVDAAMAGR   | LDGSVDFKK          | TVMVNINPEGIPVK      |
| DYWSTVK            | LEEQAQIR           | TWRNDLISATK         |
| EAQLPVIENK         | LFDSDPITVTPVEVSRK  | TYETTLK             |
| EDLIWELLNQAQEHFGK  | LGPHAGDVEGHL SFLEK | TYLPAVDEK           |
| EDTPNSVWEPK        | LGPLVEQGR          | TYLPAVDEKLR         |
| EELLPAQDIK         | LGQYASPTAK         | TYLPAVDEKLRDLYSK    |
| EESPLLIGQQSTVSDVPR | LGVRPSQGGEAPR      | TYMLAFDVNDEK        |
| EFQLFSSPHGK        | LHEAFSPVSYQHDLALLR | TYNVLD MK           |
| ELDESLQVAER        | LHIMAGR            | VDKDNEDFQESNR       |
| ELSSFIDK           | LKEEIGKELEELR      | VDVIPVNLPGEHGQR     |
| ELSYYSLEDLNNK      | LKNSLF EYQK        | VEDPESTLFGSVIR      |
| EMSGSPASGIPVK      | LKSWFEPLVEDMQR     | VFDEFKPLVEEPQNLK    |
| ENADSLQASLRPHADELK | LLDNWDSVTSTFSK     | VFSNGADLSGVTEEAPLK  |
| ENISDPTSPLR        | LLIYAVLPTGDVIGDSAK | VGDTLNLNLR          |
| EPAHLSLFGGKPMIYK   | LNAENNATFYFK       | VG FYESDVMGR        |
| EPTMYVGSTSVQTSR    | LPPNVVEESAR        | VGPEADKYR           |
| EQLGPVTQEFWDNLEK   | LQAEAFQAR          | VG YVSGWGR          |
| ESDTSYVSLK         | LQGTLPVEAR         | VIGNMGQTMEQLTPELK   |
| ESYSGVTLDPR        | LQHLENELTHDIITK    | VKDISEVVTTPR        |
| ETAVDGELVVLYDVK    | LRDLYSK            | VKDLATVYVDVLK       |
| EVAFDLEIPK         | LREQLGPVTQEFWDNLEK | VKSPELQAEAK         |
| EVDLKDYEDQQK       | LRTEGDGVYTLNDKK    | VKSPELQAEAKSYFEK    |
| EWFWDLATGTMK       | LSINTHPSQKPLSITVR  | VLEPTLK             |
| EYVLPSFEVIVEPTEK   | LSNENHGIAQR        | VLNQELR             |
| FAFNLYR            | LSPIYNLVPVK        | VLSLAQEYVGG SPEK    |
| FEDGVLPDPYPR       | LTIGEGQQHHLGGAK    | VLVDHFGYTK          |
| FFHKNEIWYR         | LVAYYTLIGASGQR     | VMDKYTFELSR         |
| FKDLGEENFK         | LVDKFLEDVKK        | VNKDDEEFIESNK       |
| FLATTPNSLLVSWQPPR  | LVNEVTEFAK         | VPEARPNM VVEHPEFLK  |
| FMETVAEK           | LVTDLTK            | VPGLYYFTYHASSR      |
| FNAVLTNPQGDYDTSTGK | LWAYLTIQELLAK      | VPGTSTSATLTGLTR     |
| FPEVDVLTK          | MAT TMIQSK         | VPLLLSEPINIIDALEM R |
| FPSGTFEQVSQLVK     | MGPTELLIEMEDWK     | VPQVSTPTLVEVSR      |
| FPVEMTHNHNFR       | MGPTELLIEMEDWKGD K | VQHIQLLQK           |
| FQNSAILTIQPK       | MKGLIDEVNQDFTNR    | VQPYLDDFQK          |
| FSVPAGIVIPSFQALTAR | MKPVPDLVPGNFK      | VRGGEGTGYFVDFSVR    |
| FSVVYAK            | MVETTAYALLTSLNLK   | VSFLSALEEYTK        |
| FSYSKNETYQLFLSYSSK | MYLG YEYVTAIR      | VTIMWTPPESAVTGYR    |
| FSYSSGHVHLSENK     | NANFKFTDHLK        | VTWAPPPSIDLTNFLVR   |
| FTNIGPDTMR         | NFPSPVDAAFR        | VVGGLVALR           |
| FTVDRPFLFLIYEHR    | NGNMAGISDQR        | VVLHPNYSQVDIGLIK    |
| FVTWIEGVMR         | NHMQYEIVIK         | VWVYPPEKK           |
| FYNQVSTPLLR        | NKPGVYTDVAYYLAWIR  | VWVYPPEKK           |
| GAYPLSIEPIGVR      | NNKDSHSLTTNIMEILR  | VYKPSAGNNSLYR       |
| GDKVWVYPPEKK       | NPANPVQR           | VYSLNDDLKPAK        |
| GDSGGAFAVQDPNDK    | NPNLPPETVDSLK      | WFYIASAFR           |
| GDSGGPLIVHK        | NSLFEYQK           | WKNFPSPVDAAFR       |

**High Confidence (continued):**

|                    |                    |                    |
|--------------------|--------------------|--------------------|
| GDSGGPLIVHKK       | NWGLSVYADKPETTK    | WLPSSSPVTGYR       |
| GDSPASSKPISINYR    | NWIQYK             | WQEEMELYR          |
| GDVAFVK            | NYNLVESLK          | WSRPQAPITGYR       |
| GETHEQVHSILHFK     | PALPAGTEDTAKEDAANR | WYEIEKIPTTFENGR    |
| GEVQAMLGQSTEELR    | PLVEEPQNLIK        | YEASILTHDSSIR      |
| GEVQAMLGQSTEELRVR  | PNSMVVEHPEFLK      | YEFLNGR            |
| GEWVALNPLR         | PPEIAHGYPEHSVR     | YEITTIHNLFR        |
| GFEPTLEALFGK       | PYTFHSHGITYYK      | YFIDFVAR           |
| GFSLDEATNLNGLLR    | QFSFPLSSEPFQGSYK   | YGLVTYATYPK        |
| GGETGYFVDFSVR      | QGHNSVFLIK         | YGMVAQVTQTLK       |
| GGETAQSadPQWEQLNNK | RHPDYSVLLLLR       | YKEENDDFASFRVDR    |
| GGYTLVSGYPK        | RHPYFYAPELLFFAK    | YLGEEYVK           |
| GHLFLQTDQPIYNPGQR  | RLDGSVDFK          | YLQEIYNSNNQK       |
| GHMLNHVER          | RLDGSVDFKK         | YLQEIYNSNNQKIVNLK  |
| GKWERPFEVK         | RLEVDIDIK          | YLYEIR             |
| GSESGIFTNTK        | RPYFPVAVGK         | YNPVVIDFEMQPIHEVLR |
| GSFEFPVGDAVSK      | RQSEDSTFYLGER      | YQISVNK            |
| GSPAINVAVHVFR      | RRDGYLFQLLR        | YTFELSR            |
| GWVTDGFSSLK        | RSFFSFLGEAFD GAR   | YVGGQEHFAHLLILR    |
| GWVTDGFSSLKDYWSTVK | RTHLPEVFLSK        | YVLPNFEVK          |
| GYSISYATK          | RVDTVDPYPYR        | YVNKEIQNAVNGVK     |
| HGGLYHENMR         | RVEPYGENFNK        | YYTYLIMNK          |
| HLEV DVWVIEPQGLR   | RYIETDPANR         |                    |

**Confident:**

|                     |                    |                    |
|---------------------|--------------------|--------------------|
| AHVDALR             | FNKPFVFLMIEQNTK    | NRDVVLTTTFVDDIK    |
| ALFVSEEEKK          | GDKVWVYPPEK        | RQDNEILIFWSK       |
| ALVEGVDQLFTDYQIK    | GETHEQVHSILHFK     | RTHLPEVFLSKVLEPTLK |
| AYYENSPQQVFSTEFVK   | GSPAINVAVHVFRK     | SDVVYTDWK          |
| DGYLFQLLR           | GSWVNKFPVEMTHNHNFR | SDVVYTDWKK         |
| DIFTGLIGPMK         | HGTDDGVVWMNWK      | SIEVFGQFNGKR       |
| DLEIEVVL FHPNYNINGK | HTSLGPLEAK         | SKEQLTPLIKK        |
| DNDGWLTSDPR         | IPKSDVVYTDWKK      | SWFEPLVEDMQR       |
| DNENVVNEYSSELEK     | IPTTFENGR          | SYFEKSKEQLTPLIKK   |
| DNEQHVK             | IQPSGGTNINEALLR    | TLLSNLEEAKK        |
| DRLDEVKEQVAEVR      | KGEWVALNPLRK       | TVIGPDGHKEVTK      |
| DSAHGFLK            | KTLLSNLEEAK        | VDTVDPYPYR         |
| DSGFQMNQLR          | LDDDLEHQGGHVLDHGHK | VQFELHYQEVK        |
| DTEEDFHVDQVTTVK     | LEQGENVFLQATDK     | VTFQLTYEEVLK       |
| DYVSQFEGSALGK       | LLPHANEVSQK        | WDPYKQGFGNVATNTDGK |
| ERGHMLNHVER         | LSPLGEEMR          | WFYIASAFRNEEYNK    |
| ESLSSYWESAK         | LSSPAVITDK         | YVGGQEHFAHLLILRDTK |
| EYHFGQAVR           | MLTPEHVFIHPGWK     |                    |

**Low Confidence:**

|               |                 |               |
|---------------|-----------------|---------------|
| DVVLTTTFVDDIK | PVWLGLGPIIK     | TEGDGVYTLNNEK |
| EDFTSLSLVLYSR | SNLDEDIIAENIVSR | VELEDWNGR     |

## INSTRUCTIONS FOR STORAGE AND USE

**Handling:** RM 8321 is a frozen aqueous solution containing approximately 440 synthetic peptides in 0.1 mL/L formic acid. Normal caution and care should be exercised during the material's handling and use.

**Storage:** The peptide mixture solution is shipped frozen (on dry ice) and, upon receipt, should be stored frozen until ready for use. A freezer temperature of  $-20\text{ }^{\circ}\text{C}$  is acceptable for storage for up to one week. If a longer storage time is anticipated, the material should be stored at or below  $-60\text{ }^{\circ}\text{C}$ . The RM should not be exposed to sunlight or ultraviolet radiation. Storage of thawed material at room or refrigerator temperatures may result in degradation or modification of constituent peptides.

**Use:** Vials of the RM to be analyzed should be removed from the freezer and allowed to stand at room temperature ( $20\text{ }^{\circ}\text{C}$  to  $25\text{ }^{\circ}\text{C}$ ) until thawed. After the material is thawed, it should be used immediately. The material should be mixed briefly with a vortex mixer before aliquots are withdrawn.

## PREPARATION AND ANALYSIS<sup>(1)</sup>

**Material Acquisition and Preparation:** The synthetic peptides used in the preparation of the RM were obtained from GenScript USA Inc. (Piscataway, NJ). Aqueous solutions of each synthetic peptide were prepared and characterized by LC-MS using a time-of-flight mass analyzer and by LC-MS/MS using an ion trap mass analyzer. Each synthetic peptide was assessed for identity, purity, and its chromatographic and mass spectrometric behavior. Based on this assessment, solutions of 440 synthetic peptides were blended together to produce the peptide mixture in RM 8321. Two additional LC-MS/MS analyses were performed on the RM to confirm the presence of each expected peptide. Using DDA, the RM was analyzed using reversed-phase LC coupled to a LTQ XL (Thermo Scientific) ion trap mass spectrometer. The data from the DDA was analyzed using theoretical fragmentation libraries of tryptic peptides from all human proteins and a library containing only the peptides used to prepare the RM. The data from the ion trap analysis was also searched using a mass spectral library containing only spectra from the peptides used to prepare the RM. In addition to the ion trap analysis, the presence of peptides in the RM was also confirmed through analysis using reversed-phase LC coupled to a model 6460 (Agilent Technologies) triple quadrupole mass spectrometer operated in dMRM mode. The dMRM method monitored three different MRM transitions for each of the expected peptides in the RM.

**Homogeneity Analysis:** The homogeneity was assessed at the time the analyses for the reference peptide content of the RM were performed. A stratified random sampling plan was devised to test for homogeneity across the production lot. The results indicated that no appreciable vial-to-vial differences were detected.

## REFERENCES

- [1] Rudnick, P.A.; Clauser, K.R.; Kilpatrick, L.E.; Tchekhovskoi, D.V.; Neta, P.; Blonder, N.; Billheimer, D.D.; Blackman, R.K.; Bunk, D.M.; Cardasis, H.L.; Ham, A.J.; Jaffe, J.D.; Kinsinger, C.R.; Mesri, M.; Neubert, T.A.; Schilling, B.; Tabb, D.L.; Tegeler, T.J.; Vega-Montoto, L.; Variyath, A.M.; Wang, M.; Wang, P.; Whiteaker, J.R.; Zimmerman, L.J.; Carr, S.A.; Fisher, S.J.; Gibson, B.W.; Paulovich, A.G.; Regnier, F.E.; Rodriguez, H.; Spiegelman, C.; Tempst, P.; Liebler, D.C.; Stein, S.E.; *Performance Metrics for Liquid Chromatography-Tandem Mass Spectrometry Systems in Proteomics Analyses*; Mol. Cell Proteomics, Vol. 9(2), pp. 225–241 (2010).
- [2] May, W.; Parris, R.; Beck II, C.; Fassett, J.; Greenberg, R.; Guenther, F.; Kramer, G.; Wise, S.; Gills, T.; Colbert, J.; Gettings, R.; MacDonald, B.; *Definition of Terms and Modes Used at NIST for Value Assignment of Reference Materials for Chemical Measurements*; NIST Special Publication 260-136 (2000); available at <http://www.nist.gov/srm/upload/SP260-136.PDF> (accessed Dec 2015).
- [3] Stein, S.; *Mass Spectral Reference Libraries: An Ever-Expanding Resource for Chemical Identification*; Anal. Chem., Vol. 84(17), pp. 7274–7282 (2012).

*Users of this RM should ensure that the Report of Investigation in their possession is current. This can be accomplished by contacting the SRM Program: telephone (301) 975-2200; fax (301) 948-3730; e-mail [srminfo@nist.gov](mailto:srminfo@nist.gov); or via the Internet at <http://www.nist.gov/srm>.*

---

<sup>(1)</sup> Certain commercial equipment, instruments, or materials are identified in this report to adequately specify the experimental procedure. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the materials or equipment identified are necessarily the best available for the purpose.